

Vaginal sample kit


This kit is for you to take samples (swabs) for testing for infections called chlamydia and gonorrhea, which are passed on from one person to another during sex.

You can have these infections without symptoms.

They are easily cured with antibiotics.

What is in my kit?

Contents may vary slightly




Your kit also includes an absorbent pad the transfer samples should be returned in and a request form that requires you to not date that you took the test.

Please refer to the STI Kit Instructions leaflet for instructions on how to return the kit to us.

How to take your vaginal swab





Remove swab from packaging without touching cotton end.



Insert cotton end of swab 2–3 inches inside vagina and gently wiggle for 15 seconds. Remove without touching the outside of your body.



Unscrew tube marked 'vaginal' without spilling liquid. Insert swab, cotton end down, and snap against lip of tube at snap-line. Screw lid back on tightly.





Find a free condoms outlet at: www.freecondoms.scot



